

## System Requirements for:

- Mediscope

with OptiPix Full, OptiPix Lite, OptiPix Dermatoscopy, OptiPix Capillaroscopy or OptiPix XYZ Measurement

Component	System Requirements
<b>Computer</b>	Intel Core i5 2410M 2,3 GHz or similar <b>Recommended systems:</b> Station: HP PRO 34XX MT, Win7P or similar Laptop: HP ProBook 4730s 17.3" HD+ LED or similar
<b>Operating System</b>	Windows XPP, Windows 7P (32 or 64 bit)
<b>Memory</b>	Minimum: 1GB <b>Recommended:</b> 2 GB
<b>Hard disk</b>	Minimum: 10 GB of Free Space <b>Recommended:</b> 100 GB of Free Space
<b>Display</b>	Laptop: 15,6 " / 1366 x 768 (WXGA) PC: 22" /1920 x 1080 (Full HD) <b>Recommended:</b> Laptop: 15,6 " / 1600x900 (HD+) PC: 24" /1920 x 1200
<b>Graphic card</b>	A Direct X compatible graphics card is recommended to avoid unnecessary load on the CPU. Please avoid using integrated simple graphics cards. Example cards are nVidia G-Force or ATI <b>Recommended:</b> ATI Radeon HD 5000/6000 Series nVIDIA GTX400/500
<b>Other</b>	Two USB2.0 ports. PCI express x1

Djupdalsvägen 22  
S-192 51, Sollentuna  
Sweden

**Phone:** +46 (0)8 35 33 60  
**Fax:** +46 (0)8 35 33 63  
**E-mail:** [Info@optilia.eu](mailto:Info@optilia.eu)  
**Web:** [www.optilia.eu](http://www.optilia.eu)

